

EDAN

**SD3**  
Ultrasonic Pocket Doppler



**SD3**  
Ultrasonic Pocket Doppler



UltraMed  
EQUIPOS MEDICOS



**EDAN** Edan Instruments, Inc.

- 3/F - B, Nanshan Medical Equipments Park, Nanhai Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
- Tel: +86-755-26898326 • Fax: +86-755-26898330 • www.edan.com.cn • Email: info@edan.com.cn
- All rights reserved. Features and specifications are subject to change without notice



ENG-OB & GYN-SD3-V1.1-20140115

Care for Health

# EDAN

## SD3

### Ultrasonic Pocket Doppler

With excellent sensitivity, complete interchangeability and high durability, SD3 series Ultrasonic Pocket Dopplers are ideal for routine fetal heart rate detection by the clinicians.



#### User-Friendly Design

- Compact design with easy one-hand operation
- Bright OLED screen with real-time FHR display
- Hi-Fi speaker

#### Interchangeable Waterproof Probe

- Plug-and-play
- 2 MHz probe with deep penetration for later pregnancy and larger patients
- 3 MHz probe detects fetal heart beat as early as 9 weeks
- 4/5/8 MHz probe for blood flow detection
- Automatic power-off while putting the probe back to the rack (SD3 PLUS & SD3 PRO)



#### Work Mode

- Real-time work mode for spot check
- Record & play with audio playback and numeric display

#### Power Source

- Standard AA batteries for easy replacements
- Built-in Li-ion rechargeable battery ensures long working time
- Unique charger stand and wall mount design



UltraMed  
EQUIPOS MEDICOS



SD3 VASCULAR



SD3 LITE



SD3



SD3 PLUS



SD3 PRO

#### Configurations

Model	SD3 VASCULAR	SD3 LITE	SD3	SD3 PLUS	SD3 PRO
FHR Display	x	x	√	√	√
Charger Stand	x	x	x	○	○
Record & Play	x	x	x	x	√
AA Batteries	○	○	○	x	x
Li-ion Battery	x	x	x	√	√

√ Standard    x Not available    ○ Optional